

Clinical Centrifuge



Brand: HWLAB

Product Code: AC012HW

Availability: 2 - 3 Days

Call for Price: +2 01000 6546

Short Description

HWLAB® Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry...


Description

HWLAB® Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry, Cellular and molecular biology, clinical Labs and blood donation centers.

features

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed form300r/min to 4000r/min
- Timer, form 15s to 99min, or continuous operation
- For 12x10ml/7ml/5ml blood collection tubes or12x15ml centrifuge tube or 8x50ml centrifuge tube
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time
- prevent sample warming and to reduce handing time



Technical specifications

Technical specification

Maximum rcf	2460xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	12x15ml/10ml/7ml/5ml tubes or 8x50ml
Speed	Form 300r/min to 4000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	370x330x175mm

Weight

11KG

Product Gallery

