

# Micropipette H model



**Brand:** HWLAB  
**Product Code:** AC008HW  
**Availability:** 2 - 3 Days  
**Call for Price:** +2 01000 6546

## Short Description

HWLAB® pipette H model is designed and manufactured to provide you with a range of robust, accurate and precise pipettes.


## Description

HWLAB® pipette H model is designed and manufactured to provide you with a range of robust, accurate and precise pipettes. The HWLAB® pipette H model is a fully adjustable, air – displacement pipette with the selected volume shown on a digital indicator (volume display).

## features

### Product features

- 3 digits
- Magnifying shape
- 121°C sterilization for low- half part
- Very low operating force
- Colour indicates pipette volume
- Positioned for perfect ergonomics
- Only a few turns to reach your desired volume



## Product Gallery

