

## ID-Card centrifuge



**Brand:** HWLAB  
**Product Code:** AC014HW  
**Availability:** 2 - 3 Days  
**Call for Price:** +2 01000 6546

### Short Description

HWLAB® ID-Card Centrifuge is widely Used in hospital,Clinical Labs...


### Description

HWLAB® ID-Card Centrifuge is widely Used in hospital,Clinical Labs and blood donation centers. This is a specially Centrifuge developed for ID-System with a capacity of 24 ID-Cards.

### features

#### Product features

- Maintenance-free brushless motor
  - LCD display of time and speed
  - Fixed speed 870r/min, rcf is 85xg
- Timer, from 15s to 99min, or continuous operation
  - Suitable for 24 ID-Cards.
- Humanized design. Three operating modes can choose ,including short spin/setting time/continuous operation
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



**Technical specifications**  
**Technical specifications**

**Maximum ref**  
**Power supply**  
**Rotor capacity**

85xg  
 AC100V-240V/50-60Hz  
 24 ID-cards

**Speed**

870r/min

**Timer**

15s-99min / continuous operation

**Motor type**

Brushless motor

**Dimension**

370x330x175mm

**Weight**

11KG

## Product Gallery

