

## Mini Clinical Centrifuge



**Brand:** HWLAB  
**Product Code:** AC011HW  
**Availability:** 2 - 3 Days  
**Call for Price:** +2 01000 6546

### Short Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis


### Description

HWLAB® mini Clinical Centrifuge is widely Used in physical and chemistry analysis, biochemistry, Cellular and molecular biology, clinical Labs and blood donation centers.

### features

#### Product features

- Maintenance-free brush less motor
- LCD display of time and speed
- Adjustable speed form 300r/min to 5000r/min
- Timer, form 15s to 99min, or continuous operation
- For 6x15ml/10ml/7ml/5ml tubes
- Lid automatically open at the end of the run to prevent sample warming and to reduce handing time



**Technical specifications**  
**Technical specifications**

<b>Maximum rcf</b>	2500xg
<b>Power supply</b>	AC100V-240V/50-60Hz
<b>Rotor capacity</b>	6x15ml/10ml/7ml/5ml tubes
<b>Speed</b>	For 300r/min to 5000r/min
<b>Timer</b>	15s-99min / continuous operation
<b>Motor type</b>	Brushless motor
<b>Dimension</b>	212x270x158mm

**Weight**

4.5KG

**Product Gallery**

