

pipet4u® pro Manual Pipettes



Brand: AHN Biotechnologie GmbH

Product Code: AC036AH

Availability: 2 - 3 Days

Call for Price: +2 01000 6546

Short Description

Reverse pipetting option for accurate and precise aspirating and dispensing of viscous liquids.

Description

Colour Coding



Ergonomic shape



Accurate and reliable pipetting

Reverse pipetting option for accurate and precise aspirating and dispensing of viscous liquids / Innovative soft grip prevents transfer of heat to internal components ensuring accurate pipette operation / Lock counter mechanism prevents volume changes

Comfortable and effortless pipetting

Light weight pipette design / Ergonomic design for easy operation / Ergonomic grip shape / Ergonomic pipetting button for easy volume setting and low strain pipetting (minimizes RSI)

Sterile handling - Fully autoclavable pipette

Easy and quick maintenance and adjustment - With the included tool

No loss of samples - Blow out system allows last drop dispensing

Compatibility



[Compatibility with AHN myTip Standard and Macro Tips](#)



[Compatibility with AHN myTip FT Filter Tips](#)



[Compatibility with AHN myTip LT Tips](#)



[Compatibility with AHN myTip LFT Filter Tips](#)

Product Gallery

